

Gentec[®] offers a complete line of Oxygen and Air Medical Gas Flowmeters to meet your clinical needs



FM187B-15L-00

- Gentec flowmeters have a 0-15 LPM flow range, with a minimum of 65 LPM flush flow rate
- The bases are available in either nickel plated solid brass or aluminum
- Durable polycarbonate hood and metering tube
- Stainless steel ball float
- Flowmeters are all calibrated at 50 psig and 70° F
- Flowmeters are all back-pressured compensated
- Color-coded for air and oxygen
- MRI (optional)
- Dual flowmeter (available)
- ISO flowmeters air/oxygen (available)
- 5-year warranty (see warranty statement for details)



FM197B-15L-DH-P

Oxygen Flowmeter w/ Power Take-off
& Chemetron Adapter

Please refer to the ordering information on the back side of the page, or contact your GENTEC[®] distributor for assistance with selecting your new GENTEC[®] flowmeters

Flowmeter Ordering Information

FLOWMETER TYPE ADAPTER	AIR FLOWMETER 0-15 LPM	AIR FLOWMETER W/ POWER TAKE-OFF	AIR DUAL FLOWMETER 0-15 LPM	OXYGEN FLOWMETER 0-15 LPM	OXYGEN W/ POWER TAKE-OFF	OXYGEN DUAL FLOWMETER 0-15 LPM
	 STANDARD INLET 1/8" NPT FEMALE					
	FM187B-15L-00	FM187B-15L-00-P	FMD187B-15L-00	FM197B-15L-00	FM197B-15L-00-P	FMD197B-15L-00
 OHMEDA ADAPTER	FM187B-15L-OH	FM187B-15L-OH-P	FMD187B-15L-OH	FM197B-15L-OH	FM197B-15L-OH-P	FMD197B-15L-OH
 CHEMETRON ADAPTER	FM187B-15L-CH	FM187B-15L-CH-P	FMD187B-15L-CH	FM197B-15L-CH	FM197B-15L-CH-P	FMD197B-15L-CH
 DISS HAND-TIGHT NUT & NIPPLE ADAPTER	FM187B-15L-DH	FM187B-15L-DH-P	FMD187B-15L-DH	FM197B-15L-DH	FM197B-15L-DH-P	FMD197B-15L-DH
 DISS NUT & NIPPLE ADAPTER	FM187B-15L-DS	FM187B-15L-DS-P	FMD187B-15L-DS	FM197B-15L-DS	FM197B-15L-DS-P	FMD197B-15L-DS
 PURITAN- BUNNETT ADAPTER	FM187B-15L-GS	FM187B-15L-GS-P	FMD187B-15L-GS	FM197B-15L-GS	FM197B-15L-GS-P	FMD197B-15L-GS
 OXEQUIP MEDSTAR ADAPTER	FM187B-15L-BS	FM187B-15L-BS-P	FMD187B-15L-BS	FM197B-15L-BS	FM197B-15L-BS-P	FMD197B-15L-BS
 OXEQUIP TWIST ADAPTER	FM187B-15L-FS	FM187B-15L-FS-P	FMD187B-15L-FS	FM197B-15L-FS	FM197B-15L-FS-P	FMD197B-15L-FS

Options:

- For aluminum flowmeters, change the B to an A (example: FM197B-15L-OH to FM197A-15L-OH)
- For dual flowmeters with 2 power take-offs, add -2P (example: FMD197B-15L-OH-2P)
- For ISO flowmeters add an E (example: FM197A-15L-OHE)
- For MRI flowmeters add MRI to the end (example: FM197A-15L-OH-MRI)

Specifications

Weight (no adapter)	Aluminum Body: 4.24 oz/122 gm Brass Body: 6.64 oz/188 gm
Increment:	0-6 LPM in ½ LPM increment, 6-15 LPM in 1 LPM Increment
Flush Flow Rate:	65 LPM (depending on supply flow)
Max Inlet Pressure:	100 psig
PTO Flow Rate:	Over 100 LPM (depending on supply flow)